


Stereoscope or Dissecting Microscope



Sep 3-3:17 PM

Microscope & Stereoscope Safety

1. Use your assigned scope.
2. Always carry your scope with 2 hands.
3. Place scope on table away from edge and watch cords.
4. When turning objectives, use 2 hands.
5. If specimens are wet or may scratch, use a paper towel on the stage. (For stereoscopes)
6. Keep objective off specimen to protect lenses.

When returning scopes,

7. Put objectives on lowest setting at lowest position.
8. Wrap cords loosely and cover.
9. Check scope back in with the teacher before putting away.

Sep 5-12:05 PM